Internet- und Körperanalyse- Waage »youw8« Internet and Body Analysis Scales

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
Receiver


Quick start guide

Use only the power supply listed in user manual

What else do you need

Computer

Wall socket

Broadband modem
or router

Warning:
Only use the cables and batteries supplied with the product.

A. Go to the Youw8.eu website

- Open a web browser on your computer.
- Go to the following url:
www.xavax.eu/youw8

B. Create a profile

- Go to Registration on this page.
- Fill in your details.
- Use an existing email address.
- Follow the instructions on the screen.

C. Enter the activation code of the white receiver

- The receiver has a sticker underneath with the activation code.
- Enter this code correctly in the required form

D. Install the receiver

- Walk to the position your modem or router is placed.
- Attach the power adapter to the receiver.
- Plug the power adapter into the wall socket.

E. Attach the receiver to your modem/router

- Plug in the Ethernet cable into the receiver.
- Plug the other end of the Ethernet cable into a free Ethernet port of your modem/router.
- The light on the receiver will first be RED, then turn YELLOW and turn GREEN after successful connection.
- If the light does not turn GREEN, please restart your modem and router and try again.

F. Select your scale and enter your activation code

- Go back to the installation screen on your PC
- Select the scale from the menu.
- The scale has also a sticker underneath with the activation code.
- Enter this activation code in the form.

G. Place the batteries in the scale

- Open the back of the scale.
- Place the batteries in the correct sequence.
- Close the compartment.
- Place the scale on a hard floor (not carpet).
- Click Next on the screen.

H. Weighing for the first time

- Follow the instructions on the screen.
- Fill-in your waist and activity level

I. Step on the scale

- Step barefoot onto the scale.
- Wait until you see your weight and then step-off.
- Go back to your PC screen.
- When you see the correct measurement in the list, select Link
- From now on the system will recognize your weights.
- On the screen you will see your Dashboard.
- You have successfully installed the product.

You now have finalized the installation and can use the product.


Caution:
Persons with medical implants (e.g. pacemakers) may not use this device, because operation of the implant could be affected.

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